
Hotplates & Stirrers  
CB300, CB302, CB500, CR300,  
CR302, SB300, SB301, SB302,  
SB500, SD300, SD500

## Instruction Manual

Version 2.3

## Before Use

**Thank you for purchasing this Stuart product. To get the best performance from the equipment, and for your own safety, please read these instructions carefully before use.**

**If the equipment is not used in the manner described in this manual the protection provided may be impaired.**

This equipment is designed to operate under the following conditions:-

- ❖ For indoor use only
- ❖ Use in a well ventilated area
- ❖ Ambient temperature range +5°C to +40°C
- ❖ Altitude to 2000m
- ❖ Relative humidity not exceeding 80%
- ❖ Mains supply fluctuation not exceeding 10%
- ❖ Overvoltage category II IEC60364-4-443
- ❖ Pollution degree 2 IEC664
- ❖ Use with a minimum distance all around of 200mm from walls or other items

## Electrical Installation

Before connection please read and understand this instruction manual and ensure that the line supply corresponds to that shown on the rating plate. All models require a supply rated at: 230V, 50Hz, ~, single phase.

Power consumption is:

| Model | Power | Model | Power |
|-------|-------|-------|-------|
| SB300 | 600W  | CB300 | 1200W |
| SB301 | 50W   | CB302 | 1250W |
| SB302 | 650W  | CB500 | 2250W |
| SB500 | 1500W | CR300 | 900W  |
| SD300 | 600W  | CR302 | 950W  |
| SD500 | 1500W |       |       |

### **THIS INSTRUMENT MUST BE EARTHED**

All models are supplied with two mains leads fitted with IEC plugs for connection to the instrument. One has a U.K. 3 pin plug and the other has a 2 pin "Schuko" plug for connection to the mains. Choose the lead appropriate for your electrical installation and discard the other.

Should neither lead be suitable, take the lead with the U.K. plug and replace the plug with a suitable alternative. This involves cutting off the moulded plug, preparing the cable and connecting to the rewirable plug in accordance with its instructions.

### **IT IS IMPORTANT THAT THIS OPERATION SHOULD ONLY BE UNDERTAKEN BY A QUALIFIED ELECTRICIAN**

NOTE: Refer to the equipment's rating plate to ensure that the plug and fusing are suitable for the voltage and wattage stated.

The wires in the mains cable are coloured as follows:

|              |           |
|--------------|-----------|
| BROWN        | - LIVE    |
| BLUE         | - NEUTRAL |
| GREEN/YELLOW | - EARTH   |

The instruments are fitted with an IEC socket at the rear of the instrument for connection of the mains lead. The appropriate mains lead should be connected BEFORE connection to the mains supply.

Should the mains lead need replacement a cable of 1mm<sup>2</sup> of harmonised code H05W-F connected to an IEC320 plug should be used.

N.B. The UK mains lead is protected by a 10A fuse mounted in the plug top.

### **IF IN DOUBT CONSULT A QUALIFIED ELECTRICIAN**

When connected to the mains electricity supply a green light on the front panel illuminates, except for the digital SD300 and SD500 models where the display will illuminate showing the word "OFF".

## **Safety Advice Before Use**

- ❖ If the mains supply is interrupted the units will restart with its current settings on the restoration of the electricity supply.
- ❖ Do not use an instrument with a ceramic top if the top plate is scratched, chipped or shows signs of chemical etching.
- ❖ Position the instrument on a firm level surface away from any heat sensitive or flammable material.
- ❖ Do not use the units with flammable liquids or in a hazardous atmosphere.
- ❖ There is a danger of liquid spills if vessels are stirred too vigorously. Always build stirrer speed slowly and never stir faster than necessary.
- ❖ Never move or carry the instrument until it has been switched off and allowed to cool.
- ❖ Never move or carry the unit with containers on the top plate or while the instrument is connected to the mains electricity supply.

## **Hot Surfaces**

**WARNING:** The top surface of the instrument may be HOT. Do not leave heaters switched on when not in use.

**WARNING:** When the surface becomes too hot to touch the red "HOT" warning light on the front panel will begin to flash. This will continue to flash while the plate temperature is above 50°C even if the unit is switched off as long as it is connected to the electricity supply.

**WARNING:** Do not disconnect or switch off the electricity supply until the "HOT" warning light has ceased flashing.

**WARNING:** The plate may still be hot when large masses are left on the plate, even though the "HOT" warning light has gone out.

**WARNING:** The top surface of the instrument may be HOT, especially in free air when a surface temperature of up to 450°C can be achieved on some of the ceramic top models and 300°C on the metal top models.

## **Heating**

### **Models SB300, SB302, SB500, CB300, CB302, CB500:**

The knob labelled "Heat" on the front panel controls the heat output and hence the plate temperature.

The knob is graduated from 1 - 9. Increasing the value increases the heat output.

Turn the hotplate control knob clockwise to the desired setting. The amber load indicator light will illuminate and it will begin to cycle on and off when the desired temperature is reached and the hotplate begins to control.

### **Models SD300 and SD500:**

When connected to the mains electricity supply the digital display will illuminate and show the word OFF.

Switch the unit on by pressing the left hand one of the 2 buttons located to the left of the digital display. The display will now show the actual plate temperature in degrees Celsius.

Press and hold the "press to set" button, the right hand one of the 2 buttons located to the left of the digital display. The display will now show the set temperature.

The set temperature may be adjusted by turning the knob located to the right of the digital display whilst continuing to hold in the "press to set" button. When the desired setting is displayed release the "press to set" button. The display will revert to show the actual plate temperature and the instrument will begin to heat.

The set temperature can be observed at any time by pressing the "press to set" button.

Turning the adjusting knob without pressing the "press to set" switch will have no effect on the hotplate setting.

#### **Models CR300 and CR302:**

The knob labelled "Heat" on the front panel controls the heat output and hence the plate temperature. When the heat is switched on the amber pilot light illuminates.

The knob is graduated from 1 - 9. Increasing the value increases the heat output.

In operation only the area defined by the pattern will become fully heated. The outer edges remain cooler for safety considerations.

Vessels larger than the described area should not be used. Failure to observe this could give rise to overheating of the outer parts of the instrument.

**NOTE: To prevent overheating of the infrared heater a temperature limiter is fitted which switches off the heater should the temperature exceed 580°C.**

## Stirring

#### **Models SB301, SB302, CB302 and CR302:**

The knob labelled "Stir" on the front panel controls the stirrer speed.

It is labelled with an arbitrary numeric scale, 1 - 9 and turning the knob to a higher setting increases stirrer speed.

When the stirrer is switched on the amber pilot light illuminates.

## Cleaning and Care

Before attempting cleaning disconnect from the mains electricity supply and ensure that the top is cool.

Periodically clean the body of the instrument using a damp cloth and mild detergent solution. Do not use harsh or abrasive cleaning agents.

**Ceramic top models:** The ceramic top is very easy to keep clean. A damp cloth will normally remove most types of contamination. For more difficult stains a domestic cream cleanser is recommended. Cleaning is made easier if spillages are attended to promptly. In any case, spillages of alkali, phosphoric acid and hydrofluoric acid **MUST** be removed immediately as these chemicals can attack and damage the glass ceramic. Ensure that the appropriate safety precautions are observed.

During cleaning and general operation take care not to scratch the surface of the top plate as this could result in subsequent thermal breakage.

**WARNING: A ceramic top which is scratched, chipped, chemically etched or damaged must not be used.**

**Metal top models:** The metal top plates should be cleaned using a damp cloth and mild detergent solution. Do not use harsh or abrasive cleaning agents. Cleaning is made easier if spillages are attended to promptly. In any case, spillages of acids and alkalis **MUST** be removed immediately as these chemicals can attack and damage the aluminium alloy. Ensure that the appropriate safety precautions are observed.

## Servicing and Repair

This product range does not require any routine servicing.

**Note: There are no internal user replaceable parts**

In the event of product failure any repair or replacement of parts should only be undertaken by suitably qualified personnel. For advice please contact the Service Department of Cole-Parmer quoting model and serial number.

Only spare parts supplied or specified by Cole-Parmer or its agents should be used. Fitting of non-approved parts may affect the performance of the safety features designed into the instrument.

**Note: The magnetic stirrer drive utilises strong magnets.**

If in any doubt, please contact the Service Department of Cole-Parmer or the point of sale.

Cole-Parmer Ltd.  
Beacon Road,  
Stone, Staffordshire,  
ST15 0SA, United Kingdom  
Tel: +44 (0)1785 810475  
Email: cpservice@coleparmer.com  
Web: www.stuart-equipment.com

## Warranty

Cole-Parmer Ltd. warrants these instruments to be free from defects in material and workmanship, when used under normal laboratory conditions, for a period of **three (3)** years. In the event of a justified claim Cole-Parmer will replace any defective component or replace the unit free of charge.

This warranty does NOT apply if damage is caused by fire, accident, misuse, neglect, incorrect adjustment or repair, damage caused by installation, adaptation, modification, fitting of non approved parts or repair by unauthorised personnel.

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**CE** This product meets the applicable EC harmonized standards for radio frequency interference and may be expected not to interfere with, or be affected by, other equipment with similar qualifications. We cannot be sure that other equipment used in its vicinity will meet these standards

and so we cannot guarantee that interference will not occur in practice. Where there is a possibility that injury, damage or loss might occur if equipment malfunctions due to radio frequency interference, or for general advice before use, contact the manufacturer.

### EU Declaration of Conformity

|                     |  |                    |      |
|---------------------|--|--------------------|------|
| <b>Product</b>      | Laboratory Equipment   | <b>File Number</b> | P225 |
| <b>Manufacturer</b> | Cole-Parmer Ltd<br>Beacon Road<br>Stone, Staffordshire<br>ST15 0SA<br>United Kingdom |                    |      |

**This declaration of conformity is issued under the sole responsibility of the manufacturer**

**Object of Declaration** Hotplates and Stirrers - Large Format  
(reference the attached list of catalogue numbers)

The object of the declaration described above is in conformity with the relevant Union Harmonisation Legislation:

|                              |            |                       |
|------------------------------|------------|-----------------------|
| <b>Low Voltage Directive</b> | 2006/95/EC | (until 19 April 2016) |
|                              | 2014/35/EU | (from 20 April 2016)  |
| <b>EMC Directive</b>         | 2014/30/EU |                       |
| <b>RoHS Directive</b>        | 2011/65/EC |                       |

References to the relevant harmonised standards used or references to the other technical specifications in relation to which conformity is declared:

|  |  |
|--|--|
| IEC/EN 61010-1:2010                            | Safety requirements for electrical equipment for measurement, control and laboratory use.<br>Part 1: General requirements.     |
| IEC/EN 61010-2-010:2014                        | Particular requirements for laboratory equipment for the heating of materials.   |
| IEC/EN 61326-1:2006<br>(SD300 and CB/SB/SD500) | Electrical equipment for measurement, control and laboratory use. EMC requirements.<br>Part 1: General requirements (Class A). |
| IEC/EN 61326-1:2013<br>(CB30X, CR30X, SB30X)   |  |

**Signed for and on behalf of the above manufacturer**

**Additional Information** Year of CE Marking: 2002

**Place of Issue** Stone, Staffordshire, UK

**Date of Issue** June 2009. Revised 16 June 2017

**Authorised Representative** Steve Marriott

**Title** Technical Director

**Signature**

# INSPECTION REPORT

MODEL \_\_\_\_\_

SERIAL NO. \_\_\_\_\_


## ELECTRICAL SAFETY

- 1. Earth continuity
- 2. Insulation
- 3. Flash test

## FUNCTIONAL

- 1. Indicators
- 2. Heater/stirrer control
- 3. Visual acceptance

QUALITY CONTROL INSPECTOR

Thank you for reading this data sheet.

For pricing or for further information, please contact us at our UK Office, using the details below.

**UK Office**

**Keison Products,**

**P.O. Box 2124, Chelmsford, Essex, CM1 3UP, England.**

**Tel: +44 (0)330 088 0560**

**Fax: +44 (0)1245 808399**

**Email: [sales@keison.co.uk](mailto:sales@keison.co.uk)**

Please note - Product designs and specifications are subject to change without notice. The user is responsible for determining the suitability of this product.