

Technical Specification

Low Temperature Disintegration Test Tumbling During Reduction



Type of furnace: Linder tester

Scope: Application – Low temperature breakdown test

of iron ore, generally in accordance with

ISO/DP 4697 (1983)

General principles: Iron Ore Pellets are placed in a rotating tube and subject to

"Tumbling During Reduction". The samples are then subjected to a sieve analysis to determine the degree of disintegration.

Retort: A heat resistant Nickel/Chrome alloy rotating reduction barrel

containing 4 lifter bars and rotating about its axis within the furnace. The dimensions of the usable section are; 130mm

internal diameter x 200mm long.

Atmosphere: A mixture of CO, CO_2 , H_2 and N_2 .

Maximum operating temperature: 1000°C (Without reduction barrel)

600°C (With reduction barrel)

Furnace zones: Single control zone.

Heated chamber dimensions: 260mm high x 375mm wide x 260mm deep

Tube axis height: 820mm

Heating elements: Wire spirals hung on ceramic support tubes mounted parallel

with the reduction barrel axis.

Insulation: Low thermal mass ceramic fibre-board pinned into a 2 piece

corrosion resistant sheet steel structure.

Construction: The insulation structure is fitted with corrosion resistant sheet

steel panels and is mounted onto a rigid rolled hollow section steel framework fitted with similar steel panels. The unit is

painted with a durable epoxy paint.

Door mechanism: The front section of the furnace is mounted on a counter

balanced parallel link mechanism, which allows easy access

for loading and unloading the reduction barrel.

Temperature sensors: 1 of type K thermocouple for control

1 off type K mineral insulated probe thermocouple for the

customer's use.

Temperature control: Eurotherm micro-processor based PID controller with digital

display of the measured value and set point.

Power control: Solid state relays

Speed control: The reduction barrel is driven from a variable speed electric

motor through a toothed rubber belt.

Location of controls: The controls are located on a panel in the main supporting

framework.

Furnace rating: 7 kw

Power supply requirements: 415volts, 50 Hz, 3 phase natural and earth, fused and isolated

at 20 amps/phase.

Gas controls, alarms & safety burn-off

Gas inlets are provided for nitrogen, carbon monoxide, carbon dioxide and hydrogen supplies in the form of bulkhead connectors on the control panel suitable for connection 10mm

pipe.

Each gas line is fitted with a flow meter, complete with control valve, a solenoid valve and non-return valve. After the non-return valve, the gas lines combine and a single pipe feeds

the reduction vessel.

Maximum gas flows are chosen to cover the requirements of

all known testing standards.

The maximum rate of flow of nitrogen is 25 litres/min.

The maximum rate of flow of carbon monoxide is 12 litres/min.

The maximum rate of flow of carbon dioxide is 10 litres/min.

The maximum rate of flow of hydrogen is 2000 cm3/min.

Approximate overall dimensions: 1500mm high x 1500mm wide x 1060mm keep (including

counter-balance weight).

Additional supply requirements: CO, CO₂, H₂ and N₂ process gas supplies of adequate purity,

isolated and regulated to nominally 5 psi and suitable for

connection to 10mm compression fittings.

CH₄ (methane) pilot burner gas isolated and regulated to

250mm wg (25mbar).

Exclusions from supply: The equipment does not include gas purification or drying

units. It is the responsibility of the end user to provide the

appropriate test gases in the appropriate qualities.



Thank you for reading this data sheet.

For pricing or for further information, please contact us at our UK Office, using the details below.

UK Office Keison Products,

P.O. Box 2124, Chelmsford, Essex, CM1 3UP, England.

Tel: +44 (0)330 088 0560 Fax: +44 (0)1245 808399

Email: sales@keison.co.uk

Please note - Product designs and specifications are subject to change without notice. The user is responsible for determining the suitability of this product.