**DB-2D with RS232
Communications
Operator's Manual**

Issue 4

Date of issue:

3rd May 2005.

Techne is a trade mark

© Techne 2002

CONTENTS

DB-2D with RS232, Operator's Manual

	<i>page</i>
SAFETY and INSTALATION	
English	4
Français	6
Deutsch	8
Español	10
THE DRI-BLOCK®	12
SPECIFICATION	13
Working Environment	14
OPERATION	15
Preparation	15
The overtemperature cutout	16
When you switch on	16
Front Panel Controls	16
RS232 Communication with the DB-2D	18
After Use	19
ADDITIONAL INFORMATION	20
Insert Blocks	20
Operator maintenance	21
Replacement parts	21
Calibration	22

NOTES

Introduction

Please read all the information in this booklet before using the unit.

Warning

HIGH TEMPERATURES ARE DANGEROUS; they can cause serious burns to operators and ignite combustible material.

Techne have taken great care in the design of these units to protect operators from hazards, but Operators should pay attention to the following points:

- USE CARE AND WEAR PROTECTIVE GLOVES TO PROTECT HANDS;
- DO NOT put hot objects on or near combustible objects;
- DO NOT operate the unit close to inflammable liquids or gases;
- DO NOT place any liquid directly in your unit;
- At all times USE COMMON SENSE.

Operator Safety

All Operators of Techne equipment must have available the relevant literature needed to ensure their safety.

It is important that only suitably trained personnel operate this equipment, in accordance with the instructions contained in this manual and with general safety standards and procedures. If the equipment is used in a manner not specified by Techne the protection provided by the equipment to the Operator may be impaired.

All Techne units have been designed to conform to international safety requirements and are fitted with an overtemperature cutout. On some models, the cutout is adjustable and should be set to suit the application. On all other models the cutout is preset to protect the unit.

If a safety problem should be encountered, switch off at the mains socket and remove the plug from the supply.

Installation

1. All Techne units are supplied with a power cable. This may be integral or plug-in.
2. Before connecting the mains supply, check the voltage against the rating plate. The rating plate is on the rear of the unit. Connect the mains cable to a suitable plug according to the table below.

Note that the unit must be earthed to ensure proper electrical safety.

<i>Connections</i>	<i>220V-240V</i>	<i>110V-120V</i>
Live	Brown	Black
Neutral	Blue	White
Earth	Green/yellow	Green

The fused plug supplied with the mains lead for use in the UK is fitted with the following value fuse to protect the cable: 5 AMP for TE-10A and TE-10D; 10AMP for the TU-20D.

The fuse in the unit protects the unit and the operator

Note that units marked 230V on the rating plate work at 220V; units marked 120V work at 110V. In both cases, however, the heating rate will degrade by approximately 8%.

3. Plug the mains cable into the socket on the rear of the unit.
4. Place the unit on a suitable bench or flat workspace, or in a fume cupboard if required, ensuring that the air inlet vents on the underside are free from obstruction.
5. Note that the following symbols may be next to the indicator lamps on the front panel of the units and have the following meanings:

~	:	the power indicator
⚡	:	the heater indicator
🔥	:	the overtemperature indicator

6. Symbols on or near the power switch of the unit have the following meanings:

I	:	mains switch On
O	:	mains switch Off

After use

When you have finished heating samples, remember that parts of the unit – the tubes, blocks and associated accessories – may be very hot. Take the precautions listed earlier.

Guarantee

The unit is guaranteed against any defect in material or workmanship for the period specified on the enclosed guarantee card. This period is from the date of purchase, and within this period all defective parts will be replaced free of charge provided that the defect is not the result of misuse, accident or negligence. Servicing under this guarantee should be obtained from the supplier.

Notwithstanding the description and specification(s) of the units contained in the Operator's Manual, Techne hereby reserves the right to make such changes as it sees fit to the units or to any component of the units.

This Manual has been prepared solely for the convenience of Techne customers and nothing in this Instruction Book shall be taken as a warranty, condition or representation concerning the description, merchantability, fitness for purpose or otherwise of the units or components.

Operator maintenance

NOTE: THAT THIS EQUIPMENT SHOULD ONLY BE DISMANTLED BY PROPERLY TRAINED PERSONNEL.

REMOVING THE SIDE, FRONT OR REAR PANELS EXPOSES POTENTIALLY LETHAL MAINS VOLTAGES.

THERE ARE NO OPERATOR MAINTAINABLE PARTS WITHIN THE EQUIPMENT.

In the unlikely event that you experience any problems with your unit which cannot easily be remedied, you should contact your supplier and return the unit if necessary. Please include any details of the fault observed and remember to return the unit in its original packing. Techne accept no responsibility for damage to units which are not properly packed for shipping: if in doubt, contact your supplier. See the Decontamination Certificate supplied with your unit.

1. Cleaning

Before cleaning your unit ALWAYS disconnect it from the power supply and allow it to cool below 50° C.

Your unit can be cleaned by wiping with a damp soapy cloth. Care should be exercised to prevent water from running inside the unit. Do not use abrasive cleaners.

2. Fuses

Your unit is protected by one or two fuses. These should only be changed by suitably qualified personnel.

If the fuses blow persistently, a serious fault is indicated and you may need to return the unit to your supplier for repair.

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THE DRI-BLOCK®

Please read all the information in this manual before using the unit.

Techne's Dri-Block® heaters provide a safe, dry, constant temperature source in the laboratory. You can use the units for incubation, boiling, inactivation, wet ashing, sample concentration, enzyme analysis and many other clinical and industrial purposes.

The units cover the temperature range from 25°C to 200°C using machined aluminium alloy blocks as the heat transfer medium: temperature control circuits are built into the unit. Independent indicator lights indicate mains power, heater action and overtemperature cut-out. An ON/OFF switch is mounted on the rear panel, together with the mains connector and a double fuse holder.

The units are constructed in PPS plastic, which is strong yet lightweight and can resist sustained temperatures as high as 260°C. Even though the units heat up rapidly, highly efficient insulation ensures that cases remains cool enough to handle even at maximum operating temperatures.


Techne's Dri-Block® range has been designed to comply with all relevant RF (Radio Frequency) interference and electrical safety regulations.

The operating instructions in this Operator's Manual covers the DB-2D with RS232 communication, part of a range of Dri-Blocks®.

The features are outlined in the 'Specification' section of this manual.

Other units covered in other manuals are:

- DB-1M • DB-2A • DB-2P • DB-2D • DB-3D
- DB-3 • DB-3A


Typical rear connections

SPECIFICATION

The DB-2D Dri-Block®

This unit has a recessed chamber which can hold one or two interchangeable insert blocks which accept test tubes or other sample containers. The Temperature is set by 3 front panel buttons with digital display.

The DB-2D is supplied with an extraction tool for removing the insert blocks from the unit and two plastic half blocks for reducing heat losses when only one aluminium block is used.

Specification

Working temperature range	Ambient +5°C to 100°C
Settable temperature range	0.0°C to 105.0°C
Temperature stability	±0.05°C at 40°C
Temperature stability	±0.1°C at 100°C
Temperature setting	3 push buttons
Temperature display	4 digit LED
Set point resolution	0.1°C
Set point accuracy	±1°C
Max temperature variation between blocks	0.2°C at 40° C
Heat up time 30° C to 37° C	8 minutes
Heat up time 30° C to 56° C	9 minutes
Heat up time 30°C to 100°C	12 minutes
Number of blocks	2

Electrical supply

<i>Voltage</i>	<i>Cycles</i>	<i>Power</i>
230V	50Hz-60Hz	400W
120V	50Hz-60Hz	400W

Working Environment

The Dri-Block[®] units are designed to work safely under the following conditions:

Ambient temp. range 5°C to 40°C

Humidity Up to 95% relative humidity, non-condensing

Note: The control specifications are quoted for an ambient temperature of 20°C for units fitted with 2 insert blocks type F3506. The specification will be closely held if the ambient temperature is in the range 10°C to 30°C. Outside this range the quoted figures may deteriorate but the unit will still work safely.


Radio frequency interference tested and passed to EN50081-1.

Immunity Tested and passed to EN50082-1

OPERATION

Preparation

1. Select the insert block(s) appropriate to the application. A list of available blocks is given on page 20 of this manual. Ensure that both the underside of the block(s) and the top of the hotplate are clean; efficient heat conduction between these two surfaces is essential.

Top view showing two insert blocks in place


2. Place the aluminium alloy insert blocks into the well of the unit and place the tubes containing the sample liquid in the blocks.
3. The heater design, temperature sensor and control circuit give good temperature control and uniformity, but make sure that there is a close fit of the tubes in the block to allow efficient heat transfer.
4. If less than the maximum number of insert blocks are fitted it is recommended that you fit plastic half blocks to fill the space either side of the aluminium block. The dummy blocks have the effect of reducing heat losses, thus improving temperature stability.
5. Plug the mains cable into the socket in the rear of the unit. Connect to the mains electricity supply with the plug provided or one wired correctly for your supply. Switch the power ON from the switch located next to the mains input lead. The mains indicator on the front panel will light (as well as the fuse neon on the mains input unit).

The Overtemperature Cut-out

In the event of a fault occurring which could cause the heater to remain on, power is removed from the heater when the block temperature reaches approximately 110°C.

When you Switch On

When you first switch on, the display will show the edition of the software which your unit uses. For example software issue "1" would be shown as follows:


It will display this for 1 second, then the actual temperature of the platten will be indicated.

The Front Panel Controls


The front panel controls consist of three buttons for controlling the display, a four digit LED display and three indicators.

The SET temperature Button


The **SET** temperature button displays the set temperature when pressed.

The UP ARROW Button

When the **SET** temperature button is held down and the **UP ARROW** button is pressed, the set temperature is increased.

DB-2D with RS232 OPERATOR'S MANUAL
The DOWN ARROW Button


When the **SET** temperature button is held down and the **DOWN ARROW** button is pressed, the set temperature is decreased.

Speed of Change of Set Temperature

Each press of the **UP ARROW** or **DOWN ARROW** buttons will increase or decrease the set temperature by 0.1°C. If the buttons are held down the temperature change will accelerate to 5° per second

Power Indicator


The top indicator shows that there is power to the unit

Heater Indicator

The next indicator shows when the heater is heating the unit. When the temperature is being set, and the new set temperature is higher than the temperature already in the unit, the heater indicator will light as the unit tries to follow the set temperature. If the light is on continuously the heater is getting constant power. The only exception is described under Over-Temperature Indicator. As the temperature approaches the set temperature the heater indicator will flash. When set temperature is reached the indicator will stay on for shorter periods. If the Dri-Block is above the set temperature then the indicator will be off, as the heater is not getting any power.

Over-Temperature Indicator

If the unit should, for any reason, exceed the temperature set for the over-temperature cutout (see page 16), the over-temperature indicator will light. The heater will have been switched off and the unit will begin to cool even if the heater light is on (the light staying on or not depends on which circuit has sensed an over-temperature).

RS232 Communication with the DB-2D

You can communicate with the DB-2D via a computer and the port on the back of the Dri-Block using the RS232 cable supplied with the Dri-Block.

Settings

The Communication port settings are 9600 baud, no parity, 8 data bits and one stop bit (9600, n, 8, 1).

Receiving temperatures

The command to receive the block temperature is the letter T and the Dri-Block will return a fixed length string: TB+1234.56<CR>LF>.

Set temperature

The command to set the block temperature is the fixed length string: S+1234.5<CR><LF>.


All transmitted characters are echoed.

Only set temperatures within the operating range of the DB-2D will be accepted.

Please note that the alpha characters are always upper case.

After Use

1. When you have finished heating samples, remember that parts of the unit – the tubes, blocks and associated accessories – may be very hot. Take the precautions listed earlier. We recommend that the blocks should be allowed to cool to 70°C before being removed from the Dri-Block®. They will still have to be handled with care
2. If you need to remove an insert block while it is hot, screw the extractor tool into the thermometer hole and lift the block out vertically. Never leave the extractor tool in the block while it is being used in the Dri-Block®.

Removing an insert block with the extractor tool

ADDITIONAL INFORMATION**Insert Blocks**

Insert blocks are made of aluminium alloy and must be ordered separately from the heater units. The following blocks are available from Techne:

<i>Part</i>	<i>Tube Size</i>	<i>Number of Holes</i>	<i>Size</i>
F3501	Plain block	None	95 x 76 x 51
F3502	6mm	30	95 x 76 x 51
F3503	10mm	20	95 x 76 x 51
F3504	12mm	20	95 x 76 x 51
F3505	13mm	20	95 x 76 x 51
F3506	15mm	12	95 x 76 x 51
F3507	16mm	12	95 x 76 x 51
F3508	19mm	8	95 x 76 x 51
F3509	25mm	6	95 x 76 x 51
F3510	10mm cuvettes	2 channels	95 x 76 x 51
F4460	Plain block	Thermometer hole only	95 x 76 x 51
F4461	9mm/7mm	10/20	95 x 76 x 51
F4462	24mm	6	95 x 76 x 51
F4463	26mm	6	95 x 76 x 51
F4464	1.5ml Eppendorf	20	95 x 76 x 51
F4465	0.5ml Eppendorf	30	95 x 76 x 51
F4466	Plastic half block	Plain	95 x 37 x 50
F4470	2 ml Eppendorf	20	95 x 76 x 51
F4471	0.2ml	72	95 x 76 x 51

Operator Maintenance

NOTE THAT THIS EQUIPMENT SHOULD ONLY BE DISMANTLED BY PROPERLY TRAINED PERSONNEL. **REMOVING THE SIDE, FRONT OR REAR PANELS EXPOSES POTENTIALLY LETHAL MAINS VOLTAGES.** THERE ARE NO OPERATOR MAINTAINABLE PARTS WITHIN THE EQUIPMENT.

In the unlikely event that you experience any problems with your Dri-Block[®] which cannot easily be remedied, you should contact your supplier and return the unit if necessary. Please include any details of the fault observed and remember to return the unit in its original packing. Techne accept no responsibility for damage to units which are not properly packed for shipping: if in doubt, contact your supplier.

1. Cleaning

Before cleaning your unit ALWAYS disconnect from the power supply and allow to cool below 50° C.

Your Dri-Block[®] can be cleaned by wiping with a damp soapy cloth. Care should be exercised to prevent water from running inside the unit. Do not use abrasive cleaners.

2. Fuses

Your unit is protected by two fuses. These should only be changed by suitably qualified personnel. If the fuses blow persistently, a serious fault is indicated and you may need to return the unit to your supplier for repair.

230V 2 x F2A; 120/100V 2 x F4A

Replacement Parts

The following parts may be obtained from your Techne dealer if replacements or alternatives are required:

<i>Part Number</i>	<i>Description</i>
6101308	Extraction tool
FCABLEUK	Mains cable with 13amp UK plug (5 amp fuse)
FCABLEEU	Mains cable with R/A Schuko plug
FCABLEUS	Mains cable with 3 pin US plug
FMW11	Techne HI-TEMP 96 [™] well plate, pack of 25
FMW12	HI-TEMP 96 [™] well cover, pack of 50
F4466	Plastic half block

Calibration

Remember that if you change the calibration from that set at the factory you may change the calibration at all temperatures. You may get different calibration with different blocks.

Make sure that you have good thermal contact between the blocks and the heater platen in the Dri-Block. In order to ensure that the calibration you are setting is correct, you will need to use an independent calibrated probe or thermometer.

Measure the actual temperature of the particular block you want to use in the Dri-Block using the calibrated probe or thermometer. If the calibration is not correct then you can follow this procedure. Make sure that the block and the unit are below 35°C.

- a Set the temperature display to 0.0°C. This will ensure that the block is at ambient temperature for the next procedure.
- b Hold down the Up and Down buttons and then press the Set button at the same time for 5 seconds. The display will change from the set temperature, 0.0, to the block temperature to indicate the start of self-calibration.
- c After a period of 17 minutes the display will change from the block temperature (eg 40.0°C) to “EEEE”. Measure the actual temperature of the block using the thermometer or probe. Press the ‘SET’ button and either the ‘UP’ or the ‘DOWN’ button to adjust the display to the same temperature as the measured value.
- d Press the ‘UP’ and ‘DOWN’ buttons together to confirm the value. The display will return to the block temperature and the heater will come full on.
- e After a period of 29.75 minutes the display will change from the block temperature (eg 95.0°C) to “EEEE”.
- f Measure the actual temperature of the block using the thermometer. Press the ‘SET’ button and either the ‘UP’ or ‘DOWN’ button to adjust the display to the same temperature as the measured value.
- g Press the ‘UP’ and the ‘DOWN’ buttons together to confirm the value. The display will return to the block temperature and the unit will control with the new calibration parameters.
- h Set the temperature to the particular temperature you require. Measure the temperature when the block temperature has stabilised. Repeat the procedure as necessary; you may need to put in an off-set.

Thank you for reading this data sheet.

For pricing or for further information, please contact us at our UK Office, using the details below.

UK Office

Keison Products,

P.O. Box 2124, Chelmsford, Essex, CM1 3UP, England.

Tel: +44 (0)330 088 0560

Fax: +44 (0)1245 808399

Email: sales@keison.co.uk

Please note - Product designs and specifications are subject to change without notice. The user is responsible for determining the suitability of this product.