Extremely compact personal vortex mixer with a low profile and small footprint for gentle mixing to vigorous re-suspension of cell or chemical pellets in 0.5 ml, 1.5 ml and 15 ml tubes.



'Continuous' or 'touch' operation

For tubes up to 20 mm diameter

Pressure sensitive cup accommodates tubes up to 20 mm diameter

Reliable and extremely quiet motor produces regulated and reproducible agitation throughout the speed range


Extremely easy to operate
- select either 'continuous' or
'touch' operation and turn the
knob to adjust the speed from
250 rpm to 3000 rpm

In 'touch' mode, agitation starts in response to pressure on the pressure sensitive cup

Exceptionally compact with a low profile and small footprint

– fits into almost any location

Rubber suction pads hold tight to the work surface and prevent the unit from 'walking' – they also absorb vibration and prevent its transmission to the workbench





Thank you for reading this data sheet.

For pricing or for further information, please contact us at our UK Office, using the details below.

UK Office Keison Products,

P.O. Box 2124, Chelmsford, Essex, CM1 3UP, England.

Tel: +44 (0)330 088 0560 Fax: +44 (0)1245 808399

Email: sales@keison.co.uk

Please note - Product designs and specifications are subject to change without notice. The user is responsible for determining the suitability of this product.