




Electrothermal Slide Dry Bench

This Slide Drying Bench aids the preparation of microscope slides at the specimen mounting stage. Accepting up to 50 slides, it has the facility for drying slides in different orientations.

Capacity

- Up to 50 slides
- The slides can be arranged on the bench in a number of ways as illustrated.

- Across drying bars
- Angled from bars
- Flat without bars



Temperature

- Up to 100°C

Pilot Light

- Clear to indicate "Power On" and Amber to indicate "Heater On"

Built-in energy regulator controlling temperature up to 100C.

Replacement elements easy fitted by user

Warranty: One year, parts and labour

Certifications: CE

Specifications		
Material	Aluminium Base and Top Stove Painted	
Max. Element Temperature	Up to 100°C	
Cat. No.	Capacity	Electrical Requirement
MH6616	50 slides	230V 50/60Hz
MH6616X1	50 slides	115v 50/60Hz
MH6616X6	50 slides	230v 50/60Hz, EU Plug



Thank you for reading this data sheet.

For pricing or for further information, please contact us at our UK Office, using the details below.



UK Office

Keison Products,

P.O. Box 2124, Chelmsford, Essex, CM1 3UP, England.

Tel: +44 (0)330 088 0560

Fax: +44 (0)1245 808399

Email: sales@keison.co.uk

Please note - Product designs and specifications are subject to change without notice. The user is responsible for determining the suitability of this product.